

Product datasheet for **HP217387**

CHURC 1 (CHURC1) Human qPCR Primer Pair (NM_145165)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	91612
Forward Sequence:	GATGGAGAAGAAATAGTTACCTATG
Reverse Sequence:	GCCGCATAACAGACACAGCATG
ACCN:	NM_145165 , NM_145165.1 , NM_145165.2 , NM_145165.3 , BC020550 , BC020550.1 , BC010517 , BM453839 , BP305870
UniProt ID:	Q8WUH1
Synonyms:	C14orf52; chch; My015
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)