

Product datasheet for **HP217389**

PIGM Human qPCR Primer Pair (NM_145167)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	93183
Forward Sequence:	CGCTTTCCTTATACCGCCTG
Reverse Sequence:	GAGGCGACAATAGAGTCCGCAT
ACCN:	<u>NM_145167</u> , <u>NM_145167.1</u> , <u>NM_145167.2</u> , <u>BC019865</u> , <u>BC019865.1</u> , <u>BC001803</u> , <u>BC011395</u> , <u>BX643162</u> , <u>NM_145167.3</u>
UniProt ID:	<u>Q9H3S5</u>
Synonyms:	GPI-MT-I
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)