

Product datasheet for **HP217431**

UBALD1 Human qPCR Primer Pair (NM_145253)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	124402
Forward Sequence:	CGCCTTTTCCAGGAGACCAAC
Reverse Sequence:	GGAGAACATGGTGAGAGCGTCA
ACCN:	BC062534 , NM_145253 , NM_145253.1 , NM_145253.2 , BC025327 , BC025327.1 , BC008681 , BC012018 , BC040160 , BU554162 , NM_145253.3
UniProt ID:	Q8TB05
Synonyms:	FAM100A; PP11303
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)