

## Product datasheet for **HP217450**

### ZNF563 Human qPCR Primer Pair (NM\_145276)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	147837
Forward Sequence:	GGAAGAACAGAATACTGAAGATCAG
Reverse Sequence:	GAGGCTAAATGTTTCTCCACACTG
ACCN:	<u><a href="#">NM_145276</a></u> , <u><a href="#">NM_145276.1</a></u> , <u><a href="#">NM_145276.2</a></u> , <u><a href="#">BC022523</a></u> , <u><a href="#">BC022523.1</a></u> , <u><a href="#">NM_145276.3</a></u>
UniProt ID:	<u><a href="#">Q8TA94</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)