

Product datasheet for **HP217477**

Rhotekin 2 (RTKN2) Human qPCR Primer Pair (NM_145307)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	219790
Forward Sequence:	GATGGAAGCCTTCTGGCAGCAT
Reverse Sequence:	CATGGTAGACTGAGGTTGCCTC
ACCN:	NM_145307 , NM_145307.1 , NM_145307.2 , NM_145307.3 , BC025765 , BC048002 , BC141821 , BC142725 , BC160167 , BG029885 , BX113509 , NM_145307.4
UniProt ID:	Q8IZC4
Synonyms:	bA531F24.1; PLEKHK1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)