

## Product datasheet for **HP217483**

### AFG1L Human qPCR Primer Pair (NM\_145315)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	246269
Forward Sequence:	AAAGCCTGTGGAACCGTTGCCG
Reverse Sequence:	CGACCTTGAGTCCTGTTGCCA
ACCN:	<u><a href="#">NM_145315</a></u> , <u><a href="#">NM_145315.3</a></u> , <u><a href="#">NM_145315.4</a></u> , <u><a href="#">BC018445</a></u> , <u><a href="#">BC018445.1</a></u> , <u><a href="#">BI754108</a></u> , <u><a href="#">BP246728</a></u> , <u><a href="#">NM_145315.5</a></u>
UniProt ID:	<u><a href="#">Q8WV93</a></u>
Synonyms:	AFG1; c222389; LACE1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)