

## Product datasheet for **HP217495**

### IL27 Human qPCR Primer Pair (NM\_145659)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	246778
Forward Sequence:	GGAATCTCACCTGCCAGGAGTG
Reverse Sequence:	TGGTGGAGATGAAGCAGAGACG
ACCN:	<u>NM_145659</u> , <u>NM_145659.1</u> , <u>NM_145659.2</u> , <u>NM_145659.3</u> , <u>BC062422</u> , <u>BC062422.1</u>
UniProt ID:	<u>Q8NEV9</u>
Synonyms:	IL-27; IL-27A; IL27A; IL27p28; IL30; p28
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)