

Product datasheet for **HP217506**

Glutathione S Transferase alpha 1 (GSTA1) Human qPCR Primer Pair (NM_145740)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2938
Forward Sequence:	GCAGACCAGAGCCATTCTCAAC
Reverse Sequence:	ACATACGGGCAGAAGGAGGATC
ACCN:	<u>NM_145740</u> , <u>NM_145740.1</u> , <u>NM_145740.2</u> , <u>NM_145740.3</u> , <u>NM_145740.4</u> , <u>BC053578</u> , <u>BC110891</u> , <u>BG430521</u> , <u>BI762344</u> , <u>BM973137</u> , <u>NM_145740.5</u>
UniProt ID:	<u>P08263</u>
Synonyms:	GST-epsilon; GST2; GSTA1-1; GTH1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)