

Product datasheet for **HP217518**

NFAM1 Human qPCR Primer Pair (NM_145912)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	150372
Forward Sequence:	TGGCACCTACTACTGCTCTGTC
Reverse Sequence:	AGCAGGAGCTTCTGTGGACTCT
ACCN:	NM_145912 , NM_145912.1 , NM_145912.2 , NM_145912.3 , NM_145912.4 , NM_145912.5 , NM_145912.6 , BC038241 , BC038241.1 , BI906855 , NM_145912.8
UniProt ID:	Q8NET5
Synonyms:	CNAIP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)