

Product datasheet for **HP217547**

NPB Human qPCR Primer Pair (NM_148896)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 256933 |
| Forward Sequence: | GGAGCTGCAACTGCACCCAG |
| Reverse Sequence: | AAGACGTTGCCTTGCCTGGT |
| ACCN: | <u>NM_148896</u> , <u>NM_148896.1</u> , <u>NM_148896.2</u> , <u>NM_148896.3</u> , <u>NM_148896.4</u> , <u>BC033173</u> , <u>BC126128</u> , <u>BC126485</u> , <u>BC146875</u> , <u>BC146890</u> , <u>BM978256</u> , <u>NM_148896.5</u> |
| UniProt ID: | <u>Q8NG41</u> |
| Synonyms: | L7; PPL7; PPNPB |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)