

Product datasheet for **HP217588**

MRFAP1L1 Human qPCR Primer Pair (NM_152301)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	114932
Forward Sequence:	GCTTATCCAGATCAAACGCAGG
Reverse Sequence:	CTCTGCCATCTTCGCAATCTCC
ACCN:	<u>NM_152301</u> , <u>NM_152301.1</u> , <u>NM_152301.2</u> , <u>NM_152301.3</u> , <u>BC008087</u> , <u>BC008087.1</u> , <u>BC066897</u> , <u>BX397546</u>
UniProt ID:	<u>Q96HT8</u>
Synonyms:	MGC9651; Morf4 family associated protein 1-like 1; PP784; PP784, MGC9651
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)