

Product datasheet for **HP217592**

C14orf172 (TRMT61A) Human qPCR Primer Pair (NM_152307)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	115708
Forward Sequence:	TCCTCTACTCCACAGACATCGC
Reverse Sequence:	CAATGGTGCGGATGATGGCGTG
ACCN:	<u>BC010167</u> , <u>NM_152307</u> , <u>NM_152307.1</u> , <u>NM_152307.2</u> , <u>BC010167.2</u> , <u>BC007241</u> , <u>BC016033</u> , <u>BQ015210</u>
UniProt ID:	<u>Q96FX7</u>
Synonyms:	C14orf172; GCD14; Gcd14p; hTRM61; TRM61
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)