

## Product datasheet for **HP217606**

### ERP27 Human qPCR Primer Pair (NM\_152321)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	121506
Forward Sequence:	GGAACACCATCTGCCTCTTTCG
Reverse Sequence:	ACTCTGTCACCATGTGGAGGCT
ACCN:	<a href="#">NM_152321</a> , <a href="#">NM_152321.1</a> , <a href="#">NM_152321.2</a> , <a href="#">NM_152321.3</a> , <a href="#">BC030218</a> , <a href="#">BC030218.1</a> , <a href="#">BE501761</a> , <a href="#">BQ574479</a> , <a href="#">NM_152321.4</a>
UniProt ID:	<a href="#">Q96DN0</a>
Synonyms:	C12orf46; PDIA8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)