

## Product datasheet for **HP217608**

### TEX29 Human qPCR Primer Pair (NM\_152324)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	121793
Forward Sequence:	ACGTGTTCTCTGCCTTGATTGTG
Reverse Sequence:	AACTCCGCCTTTTCGCTGGACT
ACCN:	<u><a href="#">NM_152324</a></u> , <u><a href="#">NM_152324.1</a></u> , <u><a href="#">NM_152324.2</a></u> , <u><a href="#">BC029889</a></u> , <u><a href="#">BC029889.1</a></u>
UniProt ID:	<u><a href="#">Q8N6K0</a></u>
Synonyms:	bA474D23.1; C13orf16
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)