

Product datasheet for **HP217614**

Acyl coenzyme A Thioesterase 4 (ACOT4) Human qPCR Primer Pair (NM_152331)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	122970
Forward Sequence:	CTTGCCACGTTGGCTCTAGCT
Reverse Sequence:	CCTAGAGAAATGCCCAAAGCCC
ACCN:	<u>BC031799</u> , <u>NM_152331</u> , <u>NM_152331.1</u> , <u>NM_152331.2</u> , <u>NM_152331.3</u> , <u>BC031799.1</u> , <u>BC090945</u> , <u>BC090945.1</u> , <u>BC117341</u> , <u>BC117343</u> , <u>BX248023</u> , <u>BX248047</u>
UniProt ID:	<u>Q8N9L9</u>
Synonyms:	PTE-Ib; PTE1B; PTE2B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)