

Product datasheet for **HP217620**

SPATA2L Human qPCR Primer Pair (NM_152339)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	124044
Forward Sequence:	TGGAGGAAGGAGTTCACCACCA
Reverse Sequence:	ACGTAGCCCATCTTCTGGAAGC
ACCN:	<u>NM_152339</u> , <u>NM_152339.1</u> , <u>NM_152339.2</u> , <u>NM_152339.3</u> , <u>BC039608</u> , <u>BC039608.1</u> , <u>BC027606</u> , <u>BC053588</u> , <u>BI917028</u> , <u>NM_152339.4</u>
UniProt ID:	<u>Q8IUW3</u>
Synonyms:	C16orf76; tamo
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)