

Product datasheet for **HP217628**

C17orf57 (EFCAB13) Human qPCR Primer Pair (NM_152347)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	124989
Forward Sequence:	CTATGATGAGCAACACGGAATGC
Reverse Sequence:	CACAGGTCAGAAACTCACCG
ACCN:	<u>NM_152347</u> , <u>NM_152347.1</u> , <u>NM_152347.2</u> , <u>NM_152347.3</u> , <u>NM_152347.4</u> , <u>BC036407</u> , <u>BC036407.2</u> , <u>BC037876</u> , <u>BC044604</u> , <u>BX471329</u> , <u>BX537519</u>
UniProt ID:	<u>Q8IY85</u>
Synonyms:	C17orf57
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)