

## Product datasheet for **HP217639**

### ZNF573 Human qPCR Primer Pair (NM\_152360)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	126231
Forward Sequence:	ATTCTCCTGCCACAGCTCCAA
Reverse Sequence:	TATGCTCTGAAGTTGTGTAGAGGA
ACCN:	<u><a href="#">BC029468</a></u> , <u><a href="#">NM_152360</a></u> , <u><a href="#">NM_152360.1</a></u> , <u><a href="#">NM_152360.2</a></u> , <u><a href="#">NM_152360.3</a></u> , <u><a href="#">BC064962</a></u> , <u><a href="#">BC064962.1</a></u> , <u><a href="#">BC042170</a></u> , <u><a href="#">BC015418</a></u> , <u><a href="#">BC051263</a></u> , <u><a href="#">BG572864</a></u> , <u><a href="#">BP274132</a></u>
UniProt ID:	<u><a href="#">Q86YE8</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)