

## Product datasheet for **HP217696**

### GPR 164 (OR51E1) Human qPCR Primer Pair (NM\_152430)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	143503
Forward Sequence:	CTGGCATGGAATCCACAGTGCT
Reverse Sequence:	TGGTGACACGAGGCAACGTAAG
ACCN:	<u><a href="#">BC022401</a>, <a href="#">NM_152430</a>, <a href="#">NM_152430.1</a>, <a href="#">NM_152430.2</a>, <a href="#">NM_152430.3</a></u>
UniProt ID:	<u><a href="#">Q8TCB6</a></u>
Synonyms:	D-GPCR; DGPCR; GPR136; GPR164; OR51E1P; OR52A3P; POGR; PSGR2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)