

## Product datasheet for **HP217708**

### FAM161B Human qPCR Primer Pair (NM\_152445)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	145483
Forward Sequence:	GCCTGGAGGAAGTGTTCAAAGC
Reverse Sequence:	CGCTGCTTCATTTCTCCAGTTC
ACCN:	<u><a href="#">NM_152445</a></u> , <u><a href="#">NM_152445.1</a></u> , <u><a href="#">NM_152445.2</a></u> , <u><a href="#">BC053909</a></u> , <u><a href="#">BC053909.1</a></u> , <u><a href="#">BI253737</a></u> , <u><a href="#">NM_152445.3</a></u>
UniProt ID:	<u><a href="#">Q96MY7</a></u>
Synonyms:	C14orf44; c14_5547
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)