

Product datasheet for **HP217745**

TLCD4 Human qPCR Primer Pair (NM_152487)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	148534
Forward Sequence:	GGAATTTTCGCCTGCTTGCAGAG
Reverse Sequence:	AGCATTGAGGCAATCCGCACGA
ACCN:	NM_152487 , NM_152487.1 , NM_152487.2 , BC019700 , BC108263 , BC137457 , BC137458 , BC139737 , BE466427 , BQ575814 , BU683936
UniProt ID:	Q96MV1
Synonyms:	TMEM56
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)