

Product datasheet for **HP217768**

ENTHD1 Human qPCR Primer Pair (NM_152512)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	150350
Forward Sequence:	ATCCTCAAGACAACCACACGGG
Reverse Sequence:	CTGTTGATGGGACAGAGTCCAG
ACCN:	<u>NM_152512, NM_152512.1, NM_152512.2, NM_152512.3, BC033895, NM_152512.4</u>
UniProt ID:	<u>Q8IYW4</u>
Synonyms:	CACNA1I; dj370M22.3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)