

Product datasheet for **HP217771**

COMMD1 Human qPCR Primer Pair (NM_152516)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	150684
Forward Sequence:	GCTGGAGAGTTGATGGCAAGTC
Reverse Sequence:	GACCTCATCAAATCCAAACACAG
ACCN:	<u>NM_152516</u> , <u>NM_152516.1</u> , <u>NM_152516.2</u> , <u>NM_152516.3</u> , <u>BC022046</u> , <u>BC022046.1</u> , <u>BC009266</u> , <u>BG182961</u> , <u>BG183486</u> , <u>BG215225</u>
UniProt ID:	<u>Q8N668</u>
Synonyms:	C2orf5; MURR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)