

Product datasheet for **HP217778**

C2CD6 Human qPCR Primer Pair (NM_152525)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	151254
Forward Sequence:	CCTCCAGAAAGAACAGATCCCG
Reverse Sequence:	AGATGGCTGGTTGGAGGATTGC
ACCN:	<u>NM_152525</u> , <u>NM_152525.1</u> , <u>NM_152525.2</u> , <u>NM_152525.3</u> , <u>NM_152525.4</u> , <u>NM_152525.5</u> , <u>BC030659</u> , <u>BC030659.2</u> , <u>BE502681</u> , <u>BG701187</u> , <u>BG773299</u> , <u>BX648732</u>
UniProt ID:	<u>Q53TS8</u>
Synonyms:	ALS2CR11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)