

Product datasheet for **HP217827**

Nesprin3 (SYNE3) Human qPCR Primer Pair (NM_152592)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	161176
Forward Sequence:	TGGCACACTGGAGACGCTACTC
Reverse Sequence:	GAGAGTGGCTTCAGGTTGTGAG
ACCN:	BC061899 , NM_152592 , NM_152592.1 , NM_152592.2 , NM_152592.3 , NM_152592.4 , BC061899.1 , BC146604 , BM701154 , BX248076 , BX342493 , BX346277 , BX648302 , NM_152592.6
UniProt ID:	Q6ZMZ3
Synonyms:	C14orf49; C14orf139; KASH3; LINC00341; NCRNA00341; Nesp3; NET53
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)