

Product datasheet for **HP217838**

ZNF540 Human qPCR Primer Pair (NM_152606)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	163255
Forward Sequence:	CCTGGACACTACCCAGAGGAAA
Reverse Sequence:	GCTCTTCCCTTGCTCCAGTAAG
ACCN:	NM_152606 , NM_152606.1 , NM_152606.2 , NM_152606.3 , NM_152606.4 , BC126478 , BC017940 , BC126480 , BE464266 , BF109653 , BX537980 , NM_152606.5
UniProt ID:	Q8NDQ6
Synonyms:	Nbla10512
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)