

Product datasheet for **HP217842**

LRRN4 Human qPCR Primer Pair (NM_152611)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	164312
Forward Sequence:	TTGGCTCCTCACGGATGCAAAG
Reverse Sequence:	AGAGTGGTGCTTTGGTCGGACT
ACCN:	<u>BC019612</u> , <u>NM_152611</u> , <u>NM_152611.1</u> , <u>NM_152611.2</u> , <u>NM_152611.3</u> , <u>NM_152611.4</u> , <u>BC027720</u> , <u>BC027720.2</u> , <u>BP252017</u> , <u>NM_152611.5</u>
UniProt ID:	<u>Q8WUT4</u>
Synonyms:	C20orf75; dj1056H1.1; NLRR-4; NLRR4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)