

Product datasheet for **HP217850**

Sphingomyelin Synthase 2 (SGMS2) Human qPCR Primer Pair (NM_152621)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	166929
Forward Sequence:	GCATTTCCAGTGTGCTCCAAAGC
Reverse Sequence:	GTAACCGTGTGACCGCTGAAGA
ACCN:	NM_152621 , NM_152621.1 , NM_152621.2 , NM_152621.3 , NM_152621.4 , NM_152621.5 , BC041369 , BC041369.2 , BC028705 , BX649142 , NM_152621.6
UniProt ID:	Q8NHU3
Synonyms:	CDL; SMS2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)