

Product datasheet for **HP217862**

CCER1 Human qPCR Primer Pair (NM_152638)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 196477 |
| Forward Sequence: | TCCAGCGTACAATAGACCGCAC |
| Reverse Sequence: | TAGAGCACACGGGTGGTTGGAA |
| ACCN: | <u>BC026075</u> , <u>NM_152638</u> , <u>NM_152638.1</u> , <u>NM_152638.2</u> , <u>NM_152638.3</u> , <u>BC026075.2</u> , <u>BC027600</u> , <u>BC027600.2</u> , <u>BC024183</u> , <u>BC034491</u> , <u>BX114695</u> , <u>NM_152638.4</u> |
| UniProt ID: | <u>Q8TC90</u> |
| Synonyms: | C12orf12 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)