

Product datasheet for **HP217917**

PATL1 Human qPCR Primer Pair (NM_152716)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	219988
Forward Sequence:	GCCAAACTGGAGCAGCCTATA
Reverse Sequence:	CACTCCGAGATGTCACAACAGC
ACCN:	<u>NM_152716</u> , <u>NM_152716.1</u> , <u>NM_152716.2</u> , <u>BC065264</u> , <u>BC065264.1</u> , <u>BC023255</u> , <u>BC048309</u> , <u>BC109038</u> , <u>BC109039</u> , <u>BC111047</u>
UniProt ID:	<u>Q86TB9</u>
Synonyms:	hPat1b; Pat1b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)