

## Product datasheet for **HP217920**

### DOK6 Human qPCR Primer Pair (NM\_152721)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	220164
Forward Sequence:	CAAGCACCTCTGCATGGAGTGT
Reverse Sequence:	CCAGGTTTGCCGTAGGCATAAG
ACCN:	<a href="#">NM_152721</a> , <a href="#">NM_152721.1</a> , <a href="#">NM_152721.2</a> , <a href="#">NM_152721.3</a> , <a href="#">NM_152721.4</a> , <a href="#">NM_152721.5</a> , <a href="#">BC096744</a> , <a href="#">BC096744.1</a> , <a href="#">BC008583</a> , <a href="#">BC019045</a> , <a href="#">BX280184</a>
UniProt ID:	<a href="#">Q6PKX4</a>
Synonyms:	DOK5L; HsT3226
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)