

Product datasheet for **HP217923**

SLC39A12 Human qPCR Primer Pair (NM_152725)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	221074
Forward Sequence:	GGCAGTATGACAGCCTCCAACA
Reverse Sequence:	GATGAGAAGGCTGCTCCTATGG
ACCN:	<u>NM_152725</u> , <u>NM_152725.1</u> , <u>NM_152725.2</u> , <u>NM_152725.3</u> , <u>BC094700</u> , <u>BC035118</u> , <u>BC047635</u> , <u>BC065917</u> , <u>BC117323</u> , <u>BC143551</u> , <u>BF963263</u>
UniProt ID:	<u>Q504Y0</u>
Synonyms:	bA570F3.1; LZT-Hs8; ZIP-12
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)