

Product datasheet for **HP217928**

TBC1D32 Human qPCR Primer Pair (NM_152730)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	221322
Forward Sequence:	TCAACAGCAGCAGGTGGCATTG
Reverse Sequence:	CCACTGGAGTAGTTCTGGGATG
ACCN:	NM_152730 , NM_152730.1 , NM_152730.2 , NM_152730.3 , NM_152730.4 , NM_152730.5 , BX112581 , NM_152730.6
UniProt ID:	Q96NH3
Synonyms:	BROMI; C6orf170; C6orf171
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)