

Product datasheet for **HP217931**

C6orf89 Human qPCR Primer Pair (NM_152734)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	221477
Forward Sequence:	AGGCTGATGTCCTTGCCCATTG
Reverse Sequence:	CTGCTCACAGTCGTTGTCCAC
ACCN:	<u>NM_152734</u> , <u>NM_152734.1</u> , <u>NM_152734.2</u> , <u>NM_152734.3</u> , <u>BC009380</u> , <u>BC021804</u> , <u>BC023605</u> , <u>BC032796</u> , <u>BC050582</u> , <u>BC095423</u> , <u>BC105946</u> , <u>BC105947</u> , <u>BC107489</u> , <u>BX538108</u> , <u>BX647220</u> , <u>NM_152734.4</u>
UniProt ID:	<u>Q6UWU4</u>
Synonyms:	BRAP; PS1TP5TP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)