

Product datasheet for **HP217948**

YTHDF3 Human qPCR Primer Pair (NM_152758)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	253943
Forward Sequence:	GCTACTTTCAAGCATACCCACCTC
Reverse Sequence:	ACAGGACATCTTCATACGGTTATTG
ACCN:	NM_152758 , NM_152758.1 , NM_152758.2 , NM_152758.3 , NM_152758.4 , NM_152758.5 , BC052970 , BC052970.1 , BC037492 , BP223018 , BP354748 , BX537996 , NM_152758.6
UniProt ID:	Q7Z739
Synonyms:	DF3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)