

Product datasheet for **HP217957**

DOS (CBARP) Human qPCR Primer Pair (NM_152769)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	255057
Forward Sequence:	TTTGAGCAGAGCCGCAAGACGC
Reverse Sequence:	CTGAGTCACACTCGTGGATGGT
ACCN:	<u>NM_152769</u> , <u>NM_152769.1</u> , <u>NM_152769.2</u> , <u>BC028156</u> , <u>BC028156.1</u> , <u>BM263341</u> , <u>BM743507</u> , <u>BX094586</u>
UniProt ID:	<u>Q8N350</u>
Synonyms:	BARP; C19orf26; DOS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)