

Product datasheet for **HP217961**

DYNLT2B Human qPCR Primer Pair (NM_152773)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	255758
Forward Sequence:	TGTGCTCAAGGAGGAAGTGGCA
Reverse Sequence:	GCACCACCATTTTGTATCGGTCA
ACCN:	<u>NM_152773</u> , <u>NM_152773.1</u> , <u>NM_152773.2</u> , <u>NM_152773.3</u> , <u>NM_152773.4</u> , <u>BC021177</u> , <u>BC015348</u> , <u>BM548106</u> , <u>BM983853</u>
UniProt ID:	<u>Q8WW35</u>
Synonyms:	SRTD17; TCTEX1D2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)