

## Product datasheet for **HP217965**

### MFSD8 Human qPCR Primer Pair (NM\_152778)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	256471
Forward Sequence:	TACAACACCAGTTTTACTTAGCGC
Reverse Sequence:	GTCTTCCTGAGTCATCCACACG
ACCN:	<u><a href="#">NM_152778</a></u> , <u><a href="#">NM_152778.1</a></u> , <u><a href="#">NM_152778.2</a></u> , <u><a href="#">BC029503</a></u> , <u><a href="#">BC029503.1</a></u> , <u><a href="#">BI553701</a></u> , <u><a href="#">BQ015234</a></u> , <u><a href="#">NM_152778.3</a></u>
UniProt ID:	<u><a href="#">Q8NHS3</a></u>
Synonyms:	CCMD; CLN7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)