

Product datasheet for **HP217972**

TAB3 Human qPCR Primer Pair (NM_152787)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	257397
Forward Sequence:	GCCCATTTTCAGTGATACCAGGC
Reverse Sequence:	GCTAACCTCTCCATCCTTGCTC
ACCN:	NM_152787 , NM_152787.1 , NM_152787.2 , NM_152787.3 , NM_152787.4 , BC032526 , BC032526.1 , BU849692 , NM_152787.5
UniProt ID:	Q8N5C8
Synonyms:	MAP3K7IP3; NAP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)