

Product datasheet for **HP217983**

COL24A1 Human qPCR Primer Pair (NM_152890)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	255631
Forward Sequence:	GAGGACTGCAAGGTGATGTTGG
Reverse Sequence:	TCCTCCAACACTGCCAGCAGTT
ACCN:	<u>NM_152890</u> , <u>NM_152890.1</u> , <u>NM_152890.2</u> , <u>NM_152890.3</u> , <u>NM_152890.4</u> , <u>NM_152890.5</u> , <u>NM_152890.6</u> , <u>BC113654</u> , <u>BC116180</u> , <u>BC117170</u> , <u>BC143715</u> , <u>BC143717</u> , <u>BC143718</u>
UniProt ID:	<u>Q17RW2</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)