

Product datasheet for **HP217991**

C22orf25 (TANGO2) Human qPCR Primer Pair (NM_152906)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	128989
Forward Sequence:	AGCACTCACCAACTACCTGCAG
Reverse Sequence:	CTTCAGGTAGGACAAGCTGTCC
ACCN:	<u>NM_152906</u> , <u>NM_152906.1</u> , <u>NM_152906.2</u> , <u>NM_152906.3</u> , <u>NM_152906.4</u> , <u>NM_152906.5</u> , <u>NM_152906.6</u> , <u>BC041339</u> , <u>BC041339.1</u> , <u>BF726547</u> , <u>BG685068</u> , <u>BI544307</u> , <u>BI911898</u> , <u>BM013784</u> , <u>BM048105</u> , <u>BP193987</u> , <u>BP303570</u> , <u>BU173566</u> , <u>BX345246</u>
UniProt ID:	<u>Q6ICL3</u>
Synonyms:	C22orf25; MECRCN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)