

## Product datasheet for **HP218006**

### OFCC1 Human qPCR Primer Pair (NM\_153003)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	266553
Forward Sequence:	GATGAGAGCAGACAAGGAGAACC
Reverse Sequence:	TCAAGGTGAGCAGTGGTGGCTT
ACCN:	<a href="#">NM_153003</a> , <a href="#">NM_153003.1</a> , <a href="#">BC113541</a> , <a href="#">BC113543</a> , <a href="#">BC143837</a>
Synonyms:	MGC142101; MGC142103; MRDS1; orofacial cleft 1 candidate 1; orofacial clefting chromosomal breakpoint region 1; OTTHUMP00000039226
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)