

## Product datasheet for **HP218020**

### Phospholipase A2 (PLB1) Human qPCR Primer Pair (NM\_153021)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	151056
Forward Sequence:	TCCAGCCGTTTCTGAGGACCTA
Reverse Sequence:	CAGCCACAACCTGGATGTCTGC
ACCN:	<u><a href="#">BC042674</a></u> , <u><a href="#">BC065041</a></u> , <u><a href="#">NM_153021</a></u> , <u><a href="#">NM_153021.1</a></u> , <u><a href="#">NM_153021.2</a></u> , <u><a href="#">NM_153021.3</a></u> , <u><a href="#">NM_153021.4</a></u> , <u><a href="#">BC042674.1</a></u> , <u><a href="#">BC065041.1</a></u> , <u><a href="#">BC150633</a></u> , <u><a href="#">BC153864</a></u> , <u><a href="#">BX109719</a></u>
UniProt ID:	<u><a href="#">Q6P1J6</a></u>
Synonyms:	PLB; PLB/LIP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)