

Product datasheet for **HP218024**

C16orf55 (SPATA33) Human qPCR Primer Pair (NM_153025)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	124045
Forward Sequence:	AGTCGGAGAAGCCTGTGGACAG
Reverse Sequence:	CGACCACTTTCTTCCTGCTGGA
ACCN:	<u>NM_153025</u> , <u>NM_153025.1</u> , <u>NM_153025.2</u> , <u>BC045670</u> , <u>BC148299</u> , <u>BC160108</u> , <u>BP345094</u> , <u>BQ575829</u> , <u>BU741532</u>
UniProt ID:	<u>Q96N06</u>
Synonyms:	C16orf55
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)