

Product datasheet for **HP218043**

TDRD9 Human qPCR Primer Pair (NM_153046)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	122402
Forward Sequence:	GTTCTGAGATGTTGCGTTGTCC
Reverse Sequence:	CGCTCAATGTCCTCAGACCAG
ACCN:	<u>NM_153046</u> , <u>NM_153046.2</u> , <u>BC016796</u> , <u>BC033662</u> , <u>BC042804</u> , <u>BC128057</u> , <u>BG719020</u> , <u>BU430972</u>
UniProt ID:	<u>Q8NDG6</u>
Synonyms:	C14orf75; HIG-1; HLS; NET54; SPGF30; SPNE
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)