

## Product datasheet for **HP218050**

### KIF19 Human qPCR Primer Pair (NM\_153209)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	124602
Forward Sequence:	TGATCGCTCACATCAGTCCTGC
Reverse Sequence:	GCGATGTGGTAGGAGACGTTC A
ACCN:	<a href="#">NM_153209</a> , <a href="#">NM_153209.1</a> , <a href="#">NM_153209.2</a> , <a href="#">NM_153209.3</a> , <a href="#">BC026362</a> , <a href="#">BC067749</a> , <a href="#">BC093631</a> , <a href="#">BC104852</a> , <a href="#">BC110989</a> , <a href="#">BC143327</a> , <a href="#">BE219924</a> , <a href="#">NM_153209.4</a>
UniProt ID:	<a href="#">Q2TAC6</a>
Synonyms:	KIF19A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)