

## Product datasheet for **HP218060**

### ZNF524 Human qPCR Primer Pair (NM\_153219)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	147807
Forward Sequence:	ACAGCATCTCGCACTCAGAGCT
Reverse Sequence:	GCATGGATGTTGCAGTGCCGC
ACCN:	<u><a href="#">NM_153219</a></u> , <u><a href="#">NM_153219.1</a></u> , <u><a href="#">NM_153219.2</a></u> , <u><a href="#">NM_153219.3</a></u> , <u><a href="#">BC067748</a></u> , <u><a href="#">BC067748.1</a></u> , <u><a href="#">BC007396</a></u> , <u><a href="#">BC014666</a></u> , <u><a href="#">BI756946</a></u> , <u><a href="#">BM790568</a></u> , <u><a href="#">NM_153219.4</a></u>
UniProt ID:	<u><a href="#">Q96C55</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)