

## Product datasheet for **HP218063**

### CEP120 Human qPCR Primer Pair (NM\_153223)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	153241
Forward Sequence:	GAGAGAAGGCATAGACTCCCAG
Reverse Sequence:	GACTGGTTGTGAGATGGAACAGG
ACCN:	<u>NM_153223</u> , <u>NM_153223.1</u> , <u>NM_153223.2</u> , <u>NM_153223.3</u> , <u>BC040527</u> , <u>BC144558</u> , <u>BX648687</u>
UniProt ID:	<u>Q8N960</u>
Synonyms:	CCDC100; JBTS31; SRTD13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)