

## Product datasheet for **HP218079**

### CEP170L Human qPCR Primer Pair (NM\_153243)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                              |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                            |
| Gene ID:          | 220766                                                                                                                                                                                                                                                                                                                                       |
| Forward Sequence: | CACCAGAAGGAAACAACGGTCG                                                                                                                                                                                                                                                                                                                       |
| Reverse Sequence: | CTGAAAGACGGATGACAGCAGC                                                                                                                                                                                                                                                                                                                       |
| ACCN:             | <u><a href="#">NM_153243</a></u>                                                                                                                                                                                                                                                                                                             |
| Synonyms:         | CEP170L; FAM68B; KIAA0470L                                                                                                                                                                                                                                                                                                                   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)                                                                                                                                                                                                                                               |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                              |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                 |



[View online »](#)