

Product datasheet for **HP218103**

DEFB119 Human qPCR Primer Pair (NM_153289)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	245932
Forward Sequence:	ATTGTCAGTCCTGCTGCCTCCA
Reverse Sequence:	CCAGCACTCAAGGTAGTCTTGG
ACCN:	<u>NM_153289</u> , <u>NM_153289.1</u> , <u>NM_153289.2</u> , <u>NM_153289.3</u> , <u>BC062212</u> , <u>BC062212.1</u> , <u>BI460707</u> , <u>BU561156</u> , <u>NM_153289.4</u>
UniProt ID:	<u>Q8N690</u>
Synonyms:	DEFB-19; DEFB-20; DEFB20; DEFB120; ESC42-RELA; ESC42-RELB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)